
HAFELE

CYLINDER CORES

- > Modular system allows the lock case to be premounted to the furniture before installation and the cylinders are added later
- > The same cylinder core will fit all the various designs of lock case
- > The combination tool on page 107 can be used to operate or test the locking mechanism of lock cores before the cylinder is added
- > Cylinder cores can be installed, repositioned or replaced quickly. Removal key is needed to remove the cylinder core
- > For correct locking direction cores turn clockwise to unlock and anti-clockwise to lock
- > Castings and cores are marked as diagram
- > When all markings are positioned identically, the doors/drawers will all operate with the same locking direction



Modular system cylinder cores are only compatible with lock cases on pages 102-105



To install: Line up lock so that key lines up with notch in cylinder housing as shown




Modular system cylinder core











Lock is now ready to use-to remove the cylinder core the removal key must be used

> 5 plate levers = 200 possible key changes	
> Zinc alloy with brass lever plates	
Accessories: 2 keys	

> Suitable for all lock cases in the Module system

Accessories: 2 keys

> Order qty: 1 pc

Finish	Nickel plated	Brass-plated
Keyed to differ	210.04.606	210.04.802
Keyed alike		
Key change 100T	210.04.615	210.04.811
Key change 101T	210.04.624	-
Key change 102T	210.04.633	-
Key change 103T	210.04.642	-
Key change 104T	210.04.651	-
Key change 105T	210.04.660	-
Key change 106T	210.04.679	-
Key change 107T	210.04.688	-
Key change 108T	210.04.697	-
Key change 109T	210.04.704	-
Master key for cylinder cores	210.02.059	-

Cylinder rosette





- > Must be ordered separately
- > Order qty: Multiples of 10 pcs

Cat. No. 219.19.970 210.04.080	Finish	Nickel-plated	Brass-plated	
Odt. 140.	Cat. No.	219.19.970	210.04.080	

Removal key

ICH TECHNOLOGY 2018, © Häfele U.K. Ltd 2018. Dimensional data not binding. We



- > For removing cylinder cores
- > Brass
- > Order qty: 1 pc

	Cat. No.
Removal key	210.02.004